# mtDNAQuant: a new molecular test for quantification of mitochondrial DNA as a read out of stress and the antioxidant effect of vitamin C

### Anna GINOSIAN<sup>1</sup>, Fabien GUILLON<sup>1</sup> and Marie-Claude AMOUREUX<sup>2,3</sup>

<sup>1</sup>Sterlab, 06224 Vallauris Cedex, FRANCE

<sup>2</sup>Eurobio, Les Ulis, 91953 Courtaboeuf Cedex B, FRANCE

<sup>3</sup>Institut des Neurosciences de la Timone UMR 7289, Aix Marseille Université, CNRS, 13385, Marseille, France

### INTRODUCTION

Molecular-based assays are useful tools to monitor cell or tissue response to a variety of target compounds, drugs, cell-signaling molecules, as well as toxic or stressful molecules. Irritating, corrosive, phototoxic, oxidative and antioxidant properties of molecules are evaluated in the field of cosmetics.

The method presented here is based on the quantification of both mitochondrial DNA (mtDNA), and genomic DNA (gDNA) by duplex quantitative PCR. The calculation of a mtDNAQuant index corresponding to the ratio mtDNA/gDNA evaluates mtDNA copy number per cell, and indirectly mitochondrial biogenesis. Cellular assays of toxicity such as hepatocytes exposed to alcohol, epidermis exposed to detergent have initially allowed us to develop this molecular tool.

In the present study, mtDNAQuant was evaluated on Sterlab reconstructed human epidermis (RHE) and human keratinocytes in culture.



# MATERIALS and METHODS

### mtQuant index is the ratio of mitochondrial DNA (VIC) and genomic DNA (Red) quantified by duplex qPCR

mtDNAQuant standard curve ————



Example of mtDNAQuant index

# **RESULTS and DISCUSSIONS**

### mtDNAQuant index of Sterlab reconstructed epidermis: Irritation & Phototoxicity tests

#### **Reconstructed human** epidermis (RHE) model

The Sterlab RHE model is a three dimensional, multilayers culture of normal human keratinocyte cells from a foreskin tissue, histologically similar to the cell layers of the human skin



#### mtDNAQuant decreases with irritating drug



mtDNAQuant decreases with UV-A+ photosensitizer

#### 12 10 Non treated Non phototoxic Photoxique (chlorpromazine chlorhydrate Histidin Photoxicity

#### mtDNAQuant decreases with UV-B+ photosensitizer



In an irritation test, mtDNA decreased by 30 % when reconstructed epidermis were exposed to a strongly irritating drug (1-bromohexane) compared to non treated reconstructed epidermis or exposed to a weakly irritating drug (Heptyl Butyrate).

Reconstructed epidermis exposed to UVB (0,3J/cm<sup>2</sup>) had a mtDNAQuant index of 42 % compared to non-treated epidermis. When exposed to UVA (6J/cm<sup>2</sup>) and a phototoxic drug (Chlorpromazine chlorhydrate), mtDNAQuant decreased by 20%.

### mtDNAQuant index on keratynocytes: oxidative stress and Vitamin C protective effect

Superoxide dismutase (SOD)

mtDNAQuant (n=6 to 9)

**Net effect of vitamin C** on mtDNAQuant (%)





Two models of oxidative stress were developed on keratinocytes in culture either by exposing cells to UVB (0,15J/cm<sup>2</sup>) or to hydrogen peroxide (H<sub>2</sub>O<sub>2</sub>). The protective effect of vitamin C was evaluated in these models using mtDNAQuant and superoxide dismutase (SOD) quantification. SOD was strongly reduced by both UV (10 fold). Vitamin C was able to induce a modest decrease of mtDNAQuant after UV exposure. H<sub>2</sub>O<sub>2</sub> lead to a 48% increase of mtDNAQuant, reflecting mitochondrial biogenesis due to stress, and vitamin C induced a significant decrease (44%) of mtDNAQuant following incubation of  $H_2O_2$ . Vitamin C was able to protect keratinocytes mitochondrial toxicity induced by  $H_2O_2$  as measured by mtDNAQuant.

# CONCLUSIONS

- mtDNAQuant index is therefore a new indicator to follow toxicity and antioxidant properties of compounds at the level of mitochondria.
- mtDNAQuant test can potentially be adapted to any cell, tissue or paradigm aimed at following mitochondrial response.

References: 1-Shokolenkko et al., 2009. Oxidative stress induces degradation of mitochondrial DNA. NucleicAcid Research, 37 (2539-2548).

2-OECD, 2013. In Vitro Skin Irritation: Reconstructed Human Epidermis Test Method. OECD Guideline for the Testing of Chemicals, N°439.

3-Degli Esposti D. et al., 2012. Mitochondrial Roles and Cytoprotection in Chronic Liver Injury. *Biochemistry Research International*, **2012**, Article ID 387626

4-Kandarova, H.et. al., 2005. ECVAM feasibility study: Can the pre-validated in vitro skin model phototoxic potency of topical phototoxins?. 5th World Congress, from http://www.mattek.com/pages/abstracts/372